



Veterinary Service & Supplies  
200 Daniels Way, Suite 240  
Freehold, NJ 07728  
Phone: 848-482-7201  
Tech Support: 908-309-6272  
Email: [veterinaryservicesandsupplies@gmail.com](mailto:veterinaryservicesandsupplies@gmail.com)  
Website: [www.vetservicesandsupplies.com](http://www.vetservicesandsupplies.com)

## VSS MOTOR PACK

**A POWERFUL UNIT WHICH MAY BE USED AFTER SCALING TEETH TO CLEAN AND POLISH. MAY BE USED FOR SECTIONING AND RESHAPING TEETH. SELECTION OF BURS ENABLES VETERINARIAN TO PERFORM ANY ENDODONTAL WORK REQUIRED.**



### THE VSS MOTOR PACK CONSISTS OF:

- A COMPACT CONTROL UNIT WHICH DELIVERS MORE POWER THAN NECESSARY FOR POLISHING TEETH
- A LIGHT-WEIGHT 5" 20,000 RPM VARIABLE SPEED MICRO MOTOR HANDPIECE
- A SNAP-ON PROPHY ANGLE
- SNAP-ON POLISHING CUPS
- CLEANER/LUBRICANT SOLUTION

### A TRAUMA KIT IS AVAILABLE AS AN ACCESSORY

The Trauma Kit consists of a latch type contra angle for holding the #2, #4, #6 & #8 round burs for cleaning out existing cavities, a #33 ½ & #35 inverted cone for shaping bur-cleaned cavities, a diamond cylinder for reshaping broken teeth, an 18mm x .45 single-sided diamond disc for cutting or sectioning, and a bur tray for cleaning burs.

**THE VSS MOTOR PACK HAS BEEN DESIGNED TO BE USED IN CONJUNCTION WITH THE  
VSS ULTRASONIC SCALER.  
3 YEAR WARRANTY**