

Like New Performance Again and Again and Again.

DuraTip® **New!**
Ultrasonic Inserts




DuraTip inserts redefine scaling by providing like new performance (with no noticeable efficiency loss) for the entire life of the insert.

In stark contrast to conventional non-DuraTip inserts which lose 25% of their scaling efficiency after only 1 mm of tip wear and an astounding 50% of their scaling efficiency after 2 mm of tip wear.

Developed over thousands of hours by engineers in close collaboration with teams of hygienists, the patented DuraTip inserts achieve the previously unthinkable: continuous excellence (consistent peak performance) through an unprecedented 3 mm of tip wear (as was validated over a year of field testing). In other words, DuraTip ultrasonic inserts are ideally suited for every stage of patient care, as their continuous like new scaling efficiency translates into lesser pinch pressures, vastly superior performance over time, and lesser chair times than is associated with partially worn conventional tips.

- **PATENTED DESIGN AND STRUCTURE:** Proprietary internal geometries work synergistically to create DuraTip's signature consistent performance and unprecedented peak scaling efficiency
- **INNOVATIVE W-STACK:** Precision-cut nickel alloy stacks are pressed into a novel W-stack configuration, which results in substantially increased handling durability and optimized vibrational energy transfer across the entirety of the insert. This translates into previously unimaginable consistency of scaling power at the tip
- **ERGONOMIC, ANTI-SLIP SOFT GRIP WITH OPTIMIZED 12 MM OUTER DIAMETER TEXTURING:** for maximum ergonomic support and ease of handling
- **100% QUALITY INSPECTION:** Each and every DuraTip insert is subjected to Parkell's multi-point quality and performance inspections by Parkell's decades-trusted quality and engineering teams, to ensure that you will have the same incredible scaling experience every time you open a DuraTip box

Scaling Efficiency after Wear (mm)*



* Standard deviation of $\pm 10\%$.

** Wear Guides for all conventional non-DuraTip inserts instruct clinicians to discard inserts once there is 2 mm of distal tip wear.

- **PERFORMANCE GUIDES:** DuraTip inserts are paired with their own color coordinated Performance Guide (included with every insert). While DuraTip inserts maintain consistent peak efficiency up to an unprecedented 3mm of wear, it is important to regularly check tip wear so that use beyond this point is minimized



SKU#	DESCRIPTION
VSS4076	VSS DuraTip® 30K Insert — Perio Slim
VSS4075	VSS DuraTip® 30K Insert — Universal Slim



VSS

For more information, contact Veterinary Services and Supplies, LLC at **848-482-7201** or E-mail veterinaryservicesandsupplies@gmail.com

Veterinary Services and Supplies, LLC